

Overbed Table

These instructions are guidelines only. Use only as instructed by your healthcare provider.

What is an Overbed Table Used For?

Overbed tables are adjustable units that are installed for comfort in reading, writing and eating in bed.

How Do I Adjust the Overbed Table?

- The overbed table should be installed at a comfortable height for the user.



WARNING Adjustable overbed tables have a powerful internal spring. Do not disassemble the table or loosen any screws. Exercise extreme caution when working near the spring.

To adjust the height (Figure 1)

- To raise the table top, place your hand underneath the table top and move it upward to the desired height. Remove hand pressure to set height.
- To lower the table top, lift the release lever and push the table top downward to the desired height. Let go of the release lever to set height.

To tilt the table top (Note: This option may not be available on your model of overbed table.)

Models with rotating knob (Figure 1)

- Twist the knob of the tilt control knob to loosen.
- Tilt the table top to the desired position. Tighten the knob to set the position.

Models with tilt release lever (Figure 2)

- Pull up on the tilt release lever.
- Tilt the table top to the desired angle and release the lever to set the position.
- Note: In the model shown, the height adjustment release lever is located on the top of the post.




Figure 1: Tilting overbed table with rotating knob

How Do I Care for the Overbed Table?

- Clean the table top with mild, soapy water.




Figure 2: Tilting overbed table with tilt release lever



Please ensure that equipment returned to the Red Cross is clean and in good condition.